

made
in
germany



Ultrasonics.Steam.Ultraclean.



Ultrasonic Cleaning Units

Elmasonic EASY

The easy way to clean and prepare samples.

- More power with Pulse function for tough contaminations and sample preparations
- Easy handling with safety shutdown for more user safety
- Versatile ultrasonic unit with high quality and value for money

Elmasonic EASY


High performance ultrasonic cleaning devices with easy operation



The 9 Elmasonic EASY cleaning units are ready for typical lab applications like sample preparation as well as for lab cleaning tasks.

Furthermore the Elmasonic EASY units are covered by the "Elma Commitment to Quality" with 2 years warranty and replacement guarantee for the first year following date of purchase.

- Activatable function Pulse is suited to remove even tenacious contaminations and for sample preparation of hardly soluble substances in the laboratory.
- Integrated Sweep function for intensive ultrasonic cleaning



Suitable accessories



The extensive accessories programme includes cooling coils, noise protection boxes and various baskets, holders and brackets.

Technical Data

| | EASY 10 / H | EASY 20 / H | EASY 30 H | EASY 40 H | EASY 60 H | EASY 100 H | EASY 120 H | EASY 180 H | EASY 300 H |
|--------------------------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|-------------|
| Tank capacity max. (l) | 0.8 | 1.75 | 2.75 | 4.25 | 5.75 | 9.5 | 12.75 | 18.0 | 28.0 |
| Tank labour capacity (l) | 0.7 | 1.2 | 1.9 | 3.2 | 4.3 | 7.5 | 11.3 | 12.9 | 20.6 |
| Ultrasonic power effective (W) | 30 | 35 | 80 | 140 | 150 | 150 | 200 | 200 | 300 |
| Heating power (W) | 60 | 120 | 200 | 200 | 400 | 400 | 800 | 800 | 1200 |
| Unit ext. dim. W/D/H (mm) | 206/120/182 | 176/176/219 | 264/176/219 | 263/176/269 | 341/196/269 | 339/281/272 | 342/285/322 | 371/347/322 | 544/347/322 |
| Tank int. dim. W/D/H (mm) | 190/85/59 | 137/123/99 | 225/124/99 | 240/118/149 | 281/132/149 | 281/222/149 | 276/216/199 | 303/273/199 | 477/274/199 |
| Basket int. dim. W/D/H (mm) | 177/73/30 | 112/103/50 | 198/106/50 | 190/105/75 | 255/115/75 | 255/200/75 | 250/190/115 | 280/250/115 | 455/250/115 |
| Weight (kg) | 2.0 | 2.1 | 3.3 | 4.0 | 5.1 | 5.9 | 7.5 | 8.5 | 11.0 |



Authorised Dealer:
RJM SYSTEMS (PTY) LTD
 sales@rjmsystems.co.za