



Pick Up **New** Products

Sanplatec **“DRY KEEPER”** (Desiccator) Series



■ Features

- This is a desiccator that removes moisture from the cabinet without maintenance.
- The auto type is always dehumidifying when the power is turned on.
- The non-auto type uses silica gel left in the cabinet for cost benefits.
- C-3BS type is a smoke type.
- The hygrometer is digital and removable.

Pick Up **New** Products

Sanplatec Desiccator Series (Aluminum Frame)



■ Features

- Material: Main body...-Acryl (anti-charged type is vinyl chloride) Frame.. Aluminum
- The autodehumidification type and nonautotype can be selected.
- Some shelves can be slid to facilitate removal of stored items.
- The spacing of the shelf plate can be arbitrarily set and the space can be maximized.
- Plate materials include acrylic, antistatic vinyl chloride plate, PET, and light-shielding type.
- It is possible to select the rubber foot type with caster.
- Doors are available in magnet type, sealed handle type, and snap lock type.