

Pick Up **New** Products

Sanplatec Vacuum Desiccator Series



■ Features

- The vacuum desiccator made of all-acryl with excellent visibility and air tightness.
- Lineup the abundant size variations.
- Optimal for storing samples anaerobic to moisture and for storing samples in an inert gas environment.
- It is used for storage of reagents and samples, degassing of liquid, storage of optical products and precision equipment, and leak testing of food packages.

■ specification

- To the -0.1Mpa with the gage pressure supplied with the pressure-resistant main unit.