

Part Number	Description	Pack	Hazmat	Hazmat2	% Carbon	% Hydrogen	% Nitrogen	% Oxygen	% Sulphur	% Vol. matter	BTU/lb	% Fixed Carbon 1 (calculated)
B2300	Coal Standard Sulphur Range 0.3%	50g	Solid						0.35			
B2301	Coal Standard Sulphur Range 0.5%	50g	Solid						0.48			
B2302	Coal Standard Sulphur Range 0.7%	50g	Solid						0.76			
B2303	Coal Standard Sulphur Range 0.89%	50g	Solid						0.87			
B2304	Coal Standard Sulphur Range 1%	50g	Solid						1.07			
B2305	Coal Standard Sulphur Range 1.5%	50g	Solid						1.46			
B2306	Coal Standard Sulphur Range 2%	50g	Solid						2.04			
B2307	Coal Standard Sulphur Range 2.5%	50g	Solid						2.63			
B2308	Coal Standard Sulphur Range 3%	50g	Solid						3.03			
B2309	Coal Standard Sulphur Range 3.5%	50g	Solid						3.54			
B2310	Coal Standard Sulphur Range 4%	50g	Solid						4.25			
B2311	Coal Standard Sulphur Range 5%	50g	Solid						5.13			
B2312	Coal Standard Sulphur Range 6%	50g	Solid						5.83			
B2313	Ultra Low Sulphur Coal Standard Sulphur Range 0%	50g	Solid						0.00			
B2314	Ultra Low Sulphur Coal Standard Sulphur Range 0.02%	50g	Solid						0.034			
B2315	Coke Standard Sulphur Range 0.61%	50g	Solid						0.61			
B2316	Coke Standard Sulphur Range 1.21%	50g	Solid						1.21			
B2317	Coke Standard Sulphur Range 0.47%	50g	Solid						0.47			
B2318	Green Petroleum Coke Standard	50g	Solid		88.46	3.11	1.45		3.66	10.97	14919	(88.62)
B2319	Green Petroleum Coke RM	50g	Solid		87.89	1.89	1.79		5.00	6.85	14204	(92.55)
B2320	Ultimate Coal CHNOS Standard (S~0.25%)	50g	Solid		69.69	5.09	1.00	(17.5)	0.41	43.54	11888	(50.15)
B2321	Ultimate Coal CHNOS Standard (S~0.75%)	50g	Solid		68.88	2.76	1.03	(5.5)	0.81	15.76	11232	(63.22)
B2322	Ultimate Coal CHNOS Standard (S~2%)	50g	Solid		45.14	3.39	0.91	(5.4)	2.27	26.3	8097	(30.81)
B2323	Ultimate Coal CHNOS Standard (S~3%)	50g	Solid		65.03	4.38	1.23	(6.99)	3.04	30.37	11571	(50.3)
B2324	Ultimate Coal CHNOS Standard (S~6%)	50g	Solid		55.70	4.10	(1.28)	(9.89)	5.39	32.24	10031	(44.12)
B2325	Ultimate Coal Standard CHNS RM Set contains 1 each of the following B2320, B2321, B2322, B2323, & B2324											
B2326	Ultra Low Sulphur Coal Standard Sulphur Range 0.05%	50g	Solid						0.055			
B2327	Ultra Low Sulphur Coal Standard Sulphur Range 0.10%	50g	Solid						0.102			
B2328	Ultra Low Sulphur Coal Standard Sulphur Range 0.15%	50g	Solid						0.146			
B2329	Ultimate Coal CHNS Standard (S~0.76%)	50g	Solid		86.20	4.36	1.20	(1.94)	0.76	17.58	14829	(76.88)
B2339	Sulphur in Coal Standard	50g	Solid						1.18			
B2340	Sulphur in Coal Standard	50g	Solid						1.71			
B2341	Sulphur in Coal Standard	50g	Solid						7.06			
B2342	Sulfur in Green Petroleum Coke	50g	Solid						0.43			
B2343	Sulfur in Green Petroleum Coke	50g	Solid						0.49			
B2344	Sulfur in Green Petroleum Coke	50g	Solid						0.88			
B2345	Sulfur in Green Petroleum Coke	50g	Solid						1.20			
B2346	Sulfur in Green Petroleum Coke	50g	Solid						2.47			
B2347	Sulfur in Green Petroleum Coke	50g	Solid						2.21			
B2348	Sulfur in Green Petroleum Coke	50g	Solid						2.50			
B2349	Sulfur in Green Petroleum Coke	50g	Solid						4.34			
B2352	Green Petroleum Coke Standard	50g	Solid		91.58	4.06	1.78		1.11	12.91	15,697	(86.95)
B2353	Green Petroleum Coke Standard	50g	Solid		95.92	1.91	0.80		0.49	5.79	14861	(94.12)
B2386	Sulfur in Green Petroleum Coke	50g	Solid						5.43			
B2395	Mercury and Chlorine in Coal	25g	Solid		79.18	5.23	(1.7)		1.50			
B2396	Mercury and Chlorine in Coal	25g	Solid						0.77			
B2397	Mercury and Chlorine in Coal	25g	Solid						0.45			
B2398	Mercury and Chlorine in Coal	25g	Solid						1.17			
B2801	Proximate Coal Standard	50g	Solid						0.41	43.54	11888	(50.15)
B2802	Proximate Coal Standard	50g	Solid						0.52	29.48	12572	(65.49)
B2803	Proximate Coal Standard	50g	Solid						0.81	15.76	11323	(63.22)
B2804	Proximate Coal Standard	50g	Solid						0.85	45.13	11278	(46.71)
B2805	Proximate Coal Standard	50g	Solid						1.42	32.27	12216	(50.47)
B2806	Proximate Coal Standard	50g	Solid						2.27	26.3	8097	(30.81)
B2809	Proximate Coal Standard	50g	Solid						3.04	30.37	11571	(50.3)
B2810	Proximate Coal Standard	50g	Solid						5.38	32.24	10031	(44.12)
B2811	Proximate Coal Standard	50g	Solid						1.80	26.39	8377	(31.97)
B2812	Proximate Coal Standard	50g	Solid						6.55	21.7	6829	(32.3)
B2813	Proximate Coal Standard	50g	Solid						0.54	7.01	11336	(72.02)
B2814	Proximate Coal Standard	50g	Solid						0.61	36.31	13182	(56.13)

Approximate values shown - certified values vary from batch to batch - refer to the analysis certificate for exact values, or contact our sales team
 Brackets denote calculated values, or values that are for reference only, not certified